

Finally... Cerebral Oximetry for *all* your patients

NeurOs Cerebral Oximetry from Mespere LifeSciences



- > Exceptional Accuracy:** Trending Accuracy of $\pm 1.5\%$
Absolute Accuracy of $\pm 4.9\%$
- > Easy to Use:** Any healthcare professional can place the sensors on the patient and obtain readings with 1 button operation
- > Affordable:** Next generation of technology with an annual service contract designed to be used for all patients.



> Specifications


StO2 measurement range..... 0 - 100%
Trending Accuracy..... $\pm 1.5\%$
Absolute Accuracy..... $\pm 4.9\%$
Resolution.....1%
Display Update Rate.....1 second
Application Site.....Forehead and Tissue

> Components

VO200 - NeurOs Sensor.....SEN-NeurOs
VO200 - NeurOs Adhesives (Single-Use).....SA-NeurOs
VO200 - NeurOs Display Software.....OS-NeurOs

> Statistical Analysis

Clinical validation of Mespere NeurOs against invasive blood sample analysis by co-oximeter



NeurOs Cerebral Oximetry is proudly manufactured by:

Mespere LifeSciences
180 Frobisher Drive
Waterloo, Ontario, N2V 2A2



Contact
Tel: (519) 884 - 7575
Fax: (519) 884 - 8810
Email: info@mespere.com